

NanoZoomer Digital Pathology NDP.view2 Access (Article 849)


To access NanoZoomer Digital Pathology Serve system as an Application:

Open the Start Menu

Click on 'NDP.view2'



Click on 'Open Location', represented by a globe icon in the top left of the application window.



A new dialogue box will come up titled 'Open Location'.

Enter location you're accessing and click 'Ok'



If you have difficulties with this or any other procedure, please call the Help Desk at 979-862-4554 or write to helpdesk@cvm.tamu.edu.

Please include the asset tag number of the computer you're working with, and a good phone number to contact you.